

Tissue Sample Shipping and Order Submission Instructions

Step 1 Create an account through our website and fill the form for your service under [Histopathology Services](#).

Step 2 Our customer service will contact you to guide you to pay online using a credit card.

Step 3 Please follow the steps to prepare your sample:

- a) Label the sample container with chemical resistant pen.
- b) The maximum size of the tissue is 20 mm x 25 mm x 4 mm.
- c) The volume of fixative should be 10-20 times of your sample (completely immersed in the solution).
- d) Put your sample in either pre-filled histology containers with 10% formalin from [Fisher](#) or in a 50 ml conical centrifuge tube filled with 10% formalin.

Step 4 Shipping:

- a) Seal the container lid with parafilm to make sure there is no leak.
- b) Put the sample container in a box with bubble wrap to make sure there is no damage during shipment.

Step 5 After the order is received, we will review your order, you will be contacted only when there is a question about the order. You will receive an email once we start to work on your sample.